



Analysis of Fullerene's using ASAP

Analysis of hydrocarbons is insensitive or even not possible using conventional ESI or APCI techniques due to the lack of proton affinity of these type molecules. In this example, Buckminsterfullerene, commonly know as C₆₀, could be easily analyzed by applying 450 °C heat from the gas stream of an APCI probe to a melting point tube containing 100 ppb of C₆₀. The ASAP technique can be used to examine other types of certain derivatized fullerenes as long as the compounds can be vaporized. Lower molecular weight more volatile impurities can also be observed.



C60_100ppb_070413105436 #23-42 RT: 0.41-0.65 AV: 20 NL: 1.94E5
 T: FTMS + p APCI corona Full ms [700.00-800.00]



Analysis of 100 ppb C60



32583 Woods Ct, Harbeson, DE 19951

Tel. (302) 945-1967

asap-ms.com